Orion Star™ Series Meter

Abridged User's Guide

- English
- Español
- Français
- Deutsch
- Italiano

ROSS and the COIL trade dress are trademarks of Thermo Fisher Scientific Inc.

AQUAfast, Cahn, Ionalyzer, ionplus, KNIpHE, No Cal, ORION, perpHect, PerpHecT, PerpHecTion, pHISA, pHuture, Pure Water, Sage, Sensing the Future, SensorLink, ROSS, ROSS Ultra, Sure-Flow, Titrator PLUS, TURBO2 and Wine Master are registered trademarks of Thermo Fisher.

1-888-pHAX-ION, A+, All in One, Aplus, AQUAsnap, AssuredAccuracy, AUTO-BAR, AUTO-CAL, AUTO DISPENSER, Auto-ID, AUTO-LOG, AUTO-READ, AUTO-STIR, Auto-Test, BOD AutoEZ, Cable-Free, CERTI-CAL, CISA, DataCOLLECT, DataPLUS, digital LogR, DirectCal, DuraProbe, Environmental Product Authority, Extra Easy/Extra Value, FAST QC, Flash Titration, Flash Titrator, GAP, GLPcal, GLPcheck, GLPdoc, ISEasy, KAP, LabConnect, LogR, Low Maintenance Triode, Minimum Stir Requirement, MSR, NISS, One-Touch, One-Touch Calibration, One-Touch Measurement, Optimum Results, Orion Star, Pentrode, pHuture MMS, pHuture Pentrode, pHuture Quatrode, pHuture Triode, Quatrode, QuiKcheK, rf link, ROSS Resolution, SAOB, SMART AVERAGING, Smart CheK, SMART STABILITY, Stacked, Star Navigator 21, Stat Face, The Enhanced Lab, ThermoSense, Triode, TRIUMpH, Unbreakable pH, Universal Access are trademarks of Thermo Fisher.

Guaranteed Success and The Technical Edge are service marks of Thermo Fisher.

PerpHecT meters are protected by U.S. patent 6,168,707.

PerpHecT ROSS are protected by U.S. patent 6,168,707.

ORION Series A meters and 900A printer are protected by U.S. patents 5,198,093, D334,208 and D346,753.

ionplus electrodes and Optimum Results solutions are protected by US Patent 5,830,338.

ROSS Ultra electrodes are protected by US patents 6,793,787.

Orion ORP Standard is protected by US Patent 6,350,367.

Orion NoCal electrodes with stabilized potential patent pending.

© Copyright 2007, Thermo Fisher Scientific Inc. All rights reserved. Question everything, and Analyze.Detect.Measure.Control are trademarks of Thermo Fisher Scientific Inc.

The specifications, descriptions, drawings, ordering information and part numbers within this document are subject to change without notice.

This publication supersedes all previous publications on this subject.

English

pg. EN-1

This User's Guide has been abridged to include the most pertinent details in English, Spanish, French, German and Italian.

*To download an extended version in your language of choice, visit our website at **www.thermo.com/water**.*

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Table of Contents

Chapter I Introduction

Password Protected Methods	EN-1
AUTO-READ™	EN-1
Stirrer Control	EN-2
SMART STABILITY™ and SMART AVERAGING™	EN-2

Chapter II Display

General Description	EN-3
---------------------------	------

Chapter III Keypad

General Description	EN-5
Convenience Features	EN-6
Audible Signals & Alarms	EN-6
Automatic Shut-off	EN-6
Icon Definitions	EN-7

Chapter IV Preparation

Installing the Power Adapter	EN-9
Installing Batteries	EN-10
Connecting the Electrodes	EN-12
Turning on the Instrument	EN-14

Chapter V Menu Setup

Navigation Tips	EN-15
General Menu Settings	EN-16
Time and Date Setup	EN-18
Continuous, Timed, or AUTO-READ™ Measurement Selection	EN-20
Selecting Measurement Mode	EN-22
Setting Calibration	EN-23

Chapter VI pH Technique

pH Setup	EN-25
pH Calibration	EN-27
pH Measurement	EN-28

Chapter VII Dissolved Oxygen Technique

DO Setup	EN-29
DO Calibration	EN-31
DO Measurement	EN-33

Chapter VIII Conductivity Technique

Conductivity Setup	EN-35
Conductivity Calibration	EN-37
Conductivity Measurement	EN-39

Chapter IX ISE Technique

ISE Setup Options	EN-41
Preparation of Standards	EN-43
ISE Calibration	EN-44
ISE Measurement	EN-45

Chapter X Terms and Conditions

Warranty	EN-47
--------------------	-------

Chapter XI Declaration of Conformity

WEEE Compliance	EN-53
---------------------------	-------

Chapter I

Introduction

Congratulations! You have selected an industry-leading Orion Star Series meter designed for electrochemistry measurement in the field or in the laboratory.

- 2 & 3-Star meters provide single measurement of pH, dissolved oxygen, or conductivity.
- 4-Star meters provide dual-parameter measurements of pH/dissolved oxygen, pH/conductivity, or pH/ISE (ion selective electrode).
- 5-Star meters provide multi-parameter measurements including mV/ORP.

Built to meet the demands of busy, multi-user laboratories or plant environments, all meters are microprocessor controlled aiding in the delivery of accurate and precise measurements. The waterproof portable meters can even withstand submersion for short periods of time without any negative effect on operation. To better meet the needs of users in environmental protection and control, food and beverage, pharmaceutical, and consumer product laboratories, the Orion Star Series meters include these key features:

- **Password Protected Methods** – Meter's memory will save up to ten custom measurements and calibrations for future reference. Password protection of each method eliminates any tampering with methods as multiple users access only the procedure most appropriate to their work.
- **AUTO-READ™** – The instrument starts a measurement and automatically prints or logs data when the reading becomes stable.

- **Stirrer Control** – Benchtop meters (3-Star and higher) have stirrer control for the 096019 Stirrer Probe and the AUTO-STIR™ BOD probe, eliminating the need for additional stir plates and stir bars.
- **SMART STABILITY™ and SMART AVERAGING™** – Remove the guesswork by automatically compensating for measurement conditions thus optimizing response time.

An easy-to-reference Quick Start Guide, attached to each meter, supports daily use.

Please read this manual thoroughly before using your benchtop or portable meter. Any use outside of these instructions could invalidate your warranty and/or cause permanent damage to the meter.


Chapter II

Display

General Description

Throughout a given process, the LCD on any Orion Star Series meter provides **Temperature** and **Calibration** data. The **setup** only appears when the meter is in setup mode. The  indicates an error condition; when displayed with the , a sensor quality issue exists.

AR will be discussed in greater detail in **Chapter V, Menu Setup**.


5 Star Meter

Depicted here is the primary LCD of the 5-Star meter capable of multi-parameter measurement.

Display

Note: the lowest 3 lines of data correspond to what is being measured. ▲

Units of measure, displayed on the right side of the screen, will flash until the reading is stable.


Chapter III

Keypad


General Description

Keypad layout is the same for all Orion Star Series meters. The portable and 2-Star benchtop meters have 9 keys. The 3-, 4-, and 5-Star benchtop meters have 10 keys due to the addition of a  button.

Portable Keypad

Benchtop Keypad

Keypad

The center area of the ergonomically designed keypad is the focal point for setting up and navigating the meter. In particular,  and  /  are used frequently to change LCD views.

Convenience features include:

- **Display Backlight** – A quick press of  will turn the backlight on and off. When the meter is operating on battery power, the backlight will automatically turn off after two minutes to conserve power. When batteries are low, the backlight will no longer turn on.
- **Automatic Shut-off** – All Orion Star Series meters will shut down after 20 minutes without a keypress. This maximizes battery power on portable meters, or benchtop meters being run on battery power.
- **Audible Signals** – The meter will beep whenever a key is pressed providing immediate verification that the user's input was received.
- **Visual Alarm Signals** – Flashing  and  indicate calibration settings need adjustment. For more detail, refer to subsequent sections of this manual discussing specific measurement techniques.

To become familiar with each function, please review these icon definitions.

Icon Definitions

Key	Description	Key	Description
	<ul style="list-style-type: none"> Turns meter on, if shut off. Toggles backlight on and off, if meter is on. If meter is on, holding down  will turn off the meter. 		<ul style="list-style-type: none"> Changes measurement mode of the selected line.
			<ul style="list-style-type: none"> Changes selected line in setup, methods and log view modes. Edits the value of the flashing digit for setup, password entry, and calibration modes.
	Switches arrow on left of the screen between 3 display lines to select and edit.		Changes selected digit to edit and moves decimal point when changing values in setup, password entry, and calibration.
	<p>Starts calibration for currently selected line in measurement mode.</p> <ul style="list-style-type: none"> If arrow is pointed to the top line and the current units are pH,  will start a pH calibration. Each time  is pressed in the calibration it will accept the current point value and move to the next cal pts until the max # of cal pts are performed, then it will return to measurement mode. 		<ul style="list-style-type: none"> Prints and logs a measurement in continuous or timed measurement modes. Prints, logs, and freezes screen when the reading becomes stable in AUTO-READ™ mode. Exits setup menu and returns to measurement mode. Accepts calibration and returns to measurement mode.
	<p>Enters setup menu starting with selected line in measurement mode:</p> <p>If arrow is pointed to top line and current units are ISE, then  will enter the ISE setup screen.</p>		Enters log view and download screens.
			Turns stirrer on and off.

NOTES

Chapter IV


Preparation

Installing Power Adapter

The universal power supply that is included with your benchtop meter is the **ONLY** power supply source recommended for use with this unit. Use of any other power supply will void your meter warranty.

The external electrical power supply is rated to be operated at 100-240 VAC, 0.5 A, 50/60 Hz.

Based on your power source, select one of the four plugs provided – 110 V, 220 V, 240V – and slide it into the grooves on the adapter. A click will be heard when the plug is properly in place.


Connect the output plug of the power supply to the power receptacle on the benchtop meter. Refer to the diagram in section **Connecting the Electrodes**.

Installing Batteries

Orion Star Series meters use four AA Alkaline batteries. Do not use lithium or rechargeable batteries. Improper installation of non-alkaline batteries could create a hazard.

Note: For benchtop meters — installation of batteries is not required if the unit will always be connected to a power source via the universal power supply. ▲


1. Confirm meter is off.
2. Gently place meter upside down on a clean, lint-free cloth to prevent scratching of LCD.
3. Remove the battery case cover.
4. Insert new batteries with the “+” pole oriented as depicted in the battery compartment housing.
5. Replace the cover.

6. Stored data, calibrations and methods will remain in the meter's nonvolatile memory when batteries are being replaced, however, date and time may need to be reset.

Note: *For portable meters* — Batteries are supplied and installed from the factory. ▲


To access the battery compartment in portable meters, loosen the two screws in the center back of the meter. Note these are captive screws that cannot be completely removed.


Connecting the Electrodes

Follow the diagrams below to correctly connect electrodes and probes to the meter. The 5-Star meter is depicted; other meter models will have fewer connections.

Benchtop Meter – Electrode Connections

Portable Meter – Electrode Connections

Some connections serve multiple uses, for example:

- Use the BNC input to connect pH, ISE and ORP combination and sensing electrodes with a waterproof BNC connector.
- Benchtop meters have a pin tip input for a separate reference electrode. These require an appropriate BNC sensing electrode for measurement.
- The 970899WP dissolved oxygen probe can be used with the BNC connector.
- Use the waterproof 8 pin MiniDIN connector for conductivity probes.
- The waterproof 8 pin MiniDIN connector is also used for automatic temperature compensation (ATC) probes.
- The DO AUTO-STIR™ probe uses the waterproof 9 pin MiniDIN connector and the stirrer jack for the smaller connector.

Turning On the Instrument

With the batteries installed in the portable meters, or either the batteries installed or line power attached to the benchtop meter, press the  to turn on the instrument. A quick press of the  after the meter has powered up will toggle the backlight on and off. When the benchtop meter is drawing line power, the backlight will stay on until turned off with the . To turn off the meter press and hold the  for 3 seconds.

Chapter V


Menu Setup

Navigation Tips

- The  enters the setup mode.

- Main menu items are shown in the top line of the LCD. To scroll through the menu, press  / . Text that will appear in the LCD is represented in the **Display** column of the tables on following pages.

- Use  to select top, middle or bottom line.
- Use  /  in combination with  to edit the values in selected line.
- Use  to accept changes and return the arrow icon to the top line.
- Use  to save changes and return to the measurement mode.

General Menu Settings

The following table describes general instrument setup.

Description	Display	Range of Values	Default Value	Method Specific
Manual Temperature Compensation Setting	gEn dEgC 25.0	-5.0 - 105	25.0	Yes
Stirrer Speed Setting (Benchtop only)	gEn Stlr 4	OFF, 1, 2, 3, 4, 5, 6, 7	4	Yes
Auto-Shut Off Selection	gEn AUtO On	OFF, On	On	No

- **Manual Temperature** setting controls temperature compensation when no temperature sensors are attached to the instrument.
- **Stirrer Speed Setting** (benchtop meters only, excluding 2-Star meter) sets the speed from 1 through 7 with speed 1 being the slowest and 7 being the fastest. A speed of 3 or 4 is recommended for most applications.
- **Auto-Shut Off** controls whether the instrument will automatically turn off the instrument after 20 minutes elapse without a key press. To disable this feature, select OFF and the instrument will continue to operate.

To access these settings:

1. In the measurement mode press
2. Press / to scroll through the setup menu until **9E n** is displayed on the top line.
3. Press to accept selection and move the arrow icon to the middle line.
4. Press / to scroll through:
 - a. **dE9C** = Manual Temperature Setting
 - b. **St Ir** = Enable and Set Stirrer Speed
 - c. **AUTO** = Auto-Shut Off enable/disable
5. Press to accept the selection and move the arrow icon to the bottom line.
 - a. Press / to identify and to edit next value.
 - b. Press to accept the selection and move the arrow icon to the top line.
6. Press to return to the measurement mode.

Time and Date Setup

The following table describes in greater detail the abbreviations shown on the display screen.

Description	Display	Range of Values	Default Value	Method Specific
Time: Hour Setting	dAtE HOUr HH12	00 - 23 Sets the hour of the current time in 24 hour format	12	No
Time: Minute Setting	dAtE tInE mm12	00 - 59 Sets the minutes for the current time	00	No
Date Format	dAtE tYPE dmY	mdY, dmY Month, day, year or day, month, year display format	mdY	No
Date: Year Setting	dAtE YEAr 2004	00 - 99 Sets the year for the current date from 2000 to 2099	04	No
Date: Month Setting	dAtE dAtE mm01	01 - 12 Month setting 01 (January) through 12 (December)	01	No
Date: Day of the Month Setting	dAtE dAY dd01	01 - 31 Day of the month	01	No

1. In the measurement mode press
2. Press / to scroll through the setup menu until *DATE* is displayed on the top line.
3. Press to confirm setting and move the arrow icon to the middle line.
4. Press / to scroll through:
 - a. *HOUR* = Current Hour
 - b. *Min* = Current Minute
 - c. *TYPE* = Format the date stamp as either "mdY" (month, day, year) or "dmY" (day, month, year)
 - d. *DATE* = Current Month (Use numeric equivalents for months of the year)
 - e. *DAY* = Current Day of the Month
 - f. *YEAR* = Current Year
5. Press to accept selection and move the arrow icon to the bottom line.
6. Press / and to edit the selected value.
7. Press to accept the selection and move the arrow icon to the top line.
8. Repeat steps 3-7 to edit time and date settings as needed.
9. Press to return to the measurement mode.

Continuous, Timed, or AUTO-READ™ Measurement Selection

The following table describes in greater detail the abbreviations shown on the display screen.

Description	Display	Range of Values	Default Value	Method Specific
Continuous, Timed, or AUTO-READ Measurement Selection	rEAd tYPE COnt	COnt, AUtO, tImE Timed, continuous, or AUTO-READ; automatically prints and logs in AUTO-READ and timed modes	AUto	Yes
Timed Reading Setting	rEAd tImE 00:00	00:05 - 99:61 Timed readings in minutes and seconds	01:00	Yes

To change settings, simply use  and  /  to move between options. Pressing  will confirm the change of a value and end the flashing display.


Orion Star Series meters allow you to select your preferred timing for measurement intervals.

- In continuous mode the instrument is constantly making measurements and updating the screen. To print and log the measurement in this mode press .
- In timed mode the instrument is constantly making measurements, updating the display, and will automatically log and print data at a pre-selected time interval.
 - In the setup menu enter the minutes and seconds between timed readings by pressing  /  to change the value of the flashing digit. Press  to change the digit to be scrolled.
 - The minimum time interval is 5 seconds. Acceptable time ranges are: **0005 - 9959**

- In AUTO-READ mode the instrument starts a measurement whenever  is pressed. Once the reading becomes stable, the display freezes while data is automatically logged and printed.
- The AUTO-READ mode will also automatically operate the **Stirrer Control**. It will stop stirring when the reading becomes stable. To set the stirrer speed, press . Use  /  to select:

9E7
5E1r
OFF-7

Selecting Measurement Mode

In the measurement mode the arrow icon on the left side of the display screen indicates which mode is selected. Use  /  to scroll through the various modes associated with a selected line. Use  to move the arrow icon to the next line, and then toggle  /  to scroll through the modes associated with the selected line.

 **pH**
mV
Rel mV
ISE
 Off

 **µS/cm** or **mS/cm** for conductivity
mg/L for TDS
ppt for Salinity
MΩ-cm for resistivity
 Off

 **%Sat** for DO percent saturation
mg/L for DO concentration
 Barometric Pressure
 Off

For setup information related to specific techniques, refer to **Chapter VI** through **Chapter IX** of this manual.

Setting Calibration

The  indicates calibration mode or calibration setup.


Before beginning a calibration, prepare the electrode or probe according to the instructions received with the electrode or probe. The calibration standards, buffers, or calibration sleeve should also be prepared prior to the calibration. After preparing the electrode or probe per manufacturer's instructions, connect the electrode or probe to the meter. Basic information on proper calibration for each measurement technique is found in the appropriate section of this manual.

General Navigation for Selecting Calibration

1. In measurement mode press  until the arrow icon on the left of the screen is pointed toward the measurement mode to be calibrated:

pH, ORP, ISE • Conductivity • DO

2. Press  /  until the appropriate icon is lit for the measurement mode to be calibrated.

3. Then press  to begin the selected calibration.

NOTES


Chapter VI

pH Technique

pH Setup

The following table describes pH setup.

Description	Display	Range of Values	Default Value	Method Specific
pH Resolution	PH rES 0.01	0.1, 0.01, 0.001	0.01	Yes
pH Buffer Set	PH bUF USA	USA (USA Buffers) 1.68, 4.01, 7.00, 10.01, 12.46 EUrO (Euro Buffers) 1.68, 4.01, 6.86, 9.18	USA	Yes

1. In the measurement mode press .
2. Press  /  to scroll through the setup menu until **PH** is displayed on the top line.
3. Press  to confirm setting and move the arrow icon to the middle line.
4. Press  /  to select **rES** for resolution.
5. Press  to select option and move the arrow icon to the bottom line.
 - a. Use  /  to set the desired resolution.

6. Press  to accept the selection and move the arrow icon to the top line.
7. Press  to move to middle line and press  /  to select **buf** for Auto-Buffer-Recognition setting and press  again to move to the bottom line.
8. Press  /  to select either **USA** or **EU-D**; press  to accept the selection.
9. Press  to return to measurement mode.

pH Calibration

1. Prepare the electrode for use according to the instructions received with the electrode.
2. In setup mode select the buffer set of NIST (*USA*) or DIN (*EUR-D*) being used for the proper buffer recognition to occur, if not set previously.
3. Press  until the arrow icon is pointing to the pH measurement line.
4. Press .
5. Rinse the electrode and ATC probe and place into the buffer.
6. Wait for the **pH** to stop flashing.
 - a. Auto-Buffer-Recognition – When the **pH** stops flashing the meter will display the temperature-corrected pH value for the buffer.
 - b. Manual Calibration – When the **pH** stops flashing the meter will display the actual buffer value read by the pH electrode. Use the  and  /  to change the pH value to the temperature-corrected pH value for the buffer.
7. Once the correct buffer value is shown on the meter display, press  to proceed to the next calibration point and repeat steps 5 through 6 or press  to save the calibration.
8. The slope will be displayed before the meter returns to the measurement mode. **SLP** is displayed in the lower field and the actual electrode slope, in percent, is displayed in the main field.
 - a. On a 1 point calibration use the  and  /  to edit the slope then press  to return to the measurement mode.
 - b. On the 2 or more point calibration the meter will automatically go to the measurement mode 2 seconds after display of **SLP**.

pH Measurement

1. Rinse the electrode in deionized water. Shake off any excess water and blot dry with lint-free tissue.
2. Place the electrode in your sample.
 - a. If you are in continuous measurement mode, the instrument will start measuring immediately. If you are using the benchtop meter and the stirrer control is disabled, then pressing  will start the stirrer. The **pH** will flash until the reading is stable. Once the reading is stable you can log and print the measurement by pressing . If the stirrer is on, press  to turn off the stirrer before removing the stirrer from the sample.
 - b. If you are in AUTO-READ™ mode then press  to start the measurement. Once the reading is stable the meter will automatically log and print the results and freeze the display. If the stirrer is enabled then the stirrer will turn on when  is pressed and turn off once the reading has stabilized.
 - c. If you are in timed measurement mode then the meter will start making measurements at the frequency you have specified during setup. It will automatically log and print each measurement. If you are using the benchtop meter and the stirrer control is enabled, then pressing  will start the stirrer. Pressing  again will turn off the stirrer.
3. Remove the electrode from the sample and rinse with deionized water, shake or blot dry, then place it in your next sample and repeat step 2.
4. Once you are finished measuring samples, rinse the electrode with deionized water and blot dry. Then consult your electrode manual for proper electrode storage instructions.

Chapter VII

Dissolved Oxygen Technique

Dissolved Oxygen Setup

The following table describes dissolved oxygen setup.

Description	Display	Range of Values	Default Value	Method Specific
Dissolved Oxygen % Saturation Resolution	dO rES 0.1	1, 0.1 DO % saturation resolution	0.1	Yes
Dissolved Oxygen Concentration Resolution	dO rES 0.01	0.1, 0.01 DO concentration "mg/L" resolution	0.01	Yes
Dissolved Oxygen Barometric Pressure Compensation Selection	dO bAr AUtO	AUto, mAn Selects internal barometer or manual pressure	AUto	Yes
Dissolved Oxygen Manual Barometric Pressure Setting	dO PrES 760.0	450.0 - 850.0 Manual pressure compensation value	760.0	Yes
Dissolved Oxygen Salinity Correction Selection	dO SAL AUto	AUto, mAn Selects method of salinity correction; meters with conductivity function only	AUto	Yes
Dissolved Oxygen Manual Salinity Correction Factor	dO SALF 0	0 - 45 Manual salinity correction factor	0	Yes
Dissolved Oxygen Calibration Type Selection	dO CALt Alr	Alr, H2O, mAn, SEt0 Alr = Water sat. air H2O = Air sat. water mAn = Manual SEt0 = Zero point	Alr	Yes

D0 Setup Steps

1. From the measurement mode, press .
2. Press  /  to scroll through the setup menu until **dO** is displayed on the top line.
3. Press  to confirm the setting and move the arrow icon to the middle line.
4. Press  /  to scroll through:
 - a. **rES + % Sat** = % Saturation Resolution
 - b. **rES + mg/L** = Concentration Resolution
 - c. **bAr** = Barometer Type (Auto/Manual)
 - d. **P-rES** = Manual Pressure Compensation
 - e. **SAL** = Automatic/Manual Salinity Compensation
 - f. **SALF** = Manual Salinity Correction
 - g. **CALt** = Calibration Type
5. Press  to select the option and move the arrow icon to the bottom line.
6. Press  /  and  to enter the value. The previous table identifies what you will see on the display and the range of values which may be entered.
7. After entering a value, press  to complete programming of that option and move the arrow icon to the top line. Repeat steps 3 through 6 for all options.
8. Press  to return to measurement mode.

DO Calibration

Prior to calibration, the probe must be prepared and polarized.

- The DO probe is continuously polarized when connected to the meter. When first connected, or if more than 60 minutes has elapsed with the probe disconnected, re-connect the probe and allow 30 to 60 minutes for polarization. If the probe readings are stable, probe disconnections of less than one hour will require 5 to 25 minutes for polarization.

Note: The Star Series meters will supply a polarization current to the DO probe even with the power off. To maximize battery life in the meter, unplug the DO probe if it will not be used for an extended period. ▲

- Zeroing the probe – a DO probe can generate a 0.02 to 0.05 mg/L positive error in an oxygen free (anoxic) solution. If this error is unacceptable, then zero the probe when using a new sensing membrane, using fresh filling solution, or when measuring dissolved oxygen below 1 mg/L or 10% saturation.

An air calibration should be done prior to the zero calibration.

- In the measurement mode press .
- Press  /  until **dO** is shown on the top line.
- Press  to confirm the setting and move the arrow to the middle line.
- Press  /  to scroll to **CALC**.
- Press  to select setting and move arrow icon to the bottom line.
- Press  /  to select one of the following calibration modes.
 - Air** denotes calibration performed in water saturated air using the air calibration sleeve. This is the simplest and most accurate calibration method.

It is also the meter's default setting. If the calibration selection in  has not been changed, then pressing  will automatically perform an air calibration. Due to the inherent differences between water-saturated air and air-saturated water, upon stability, the air standard is set to 102.3% saturation.

- i. The highest possible accuracy is reached when calibration temperature is the same as the measuring temperature.
 - ii. Moisten the sponge or absorbent cloth in the calibration sleeve with distilled water and insert the probe into the sleeve, but without touching the water saturated material. For BOD measurements, this calibration can be performed in a BOD bottle.
 - iii. For oxygen levels below 1 mg/L, a zero calibration point is often required.
- b. **H₂O** denotes calibration performed using a water sample that is 100% saturated with air. It is the least commonly used calibration method.
- c. **mAn** denotes manual calibration using a water sample with a known concentration of dissolved oxygen. It can be used to calibrate the sensor to the value achieved by a Winkler titration.
- i. A Winkler calibration involves performing a manual Winkler titration and then using that sample as a standard. The titration oxygen level result is entered in a Winkler calibration as the DO standard value. This correlates the meter input to the Winkler titration. Note, this method is inherently less accurate due to the possibility of titration errors introduced when the calibration is set to the titration test results.
- d. **SEtO** denotes a zero calibration, which is used for very low level DO measurements. It is not generally required unless you are making measurements below 5% Saturation or 0.5 mg/L.

7. Press  to accept selection and return the arrow icon to the top line.
8. Press  to return to measurement mode.
9. The probe and calibration standard (water-saturated air or air-saturated water) should be allowed to reach equilibrium before the system is calibrated.
 - a. Press .
 - b. Wait for reading to stabilize.
 - c. Meter will display 102.3% and return to measurement mode.

DO Measurement

1. Rinse the DO probe, as well as stirrer if being used, in deionized water. Blot dry with a lint-free tissue.
2. Place the DO probe in your sample.
 - a. If you are in continuous measurement mode, the instrument will start measuring immediately. If you are using the benchtop meter and the stirrer control is enabled, then pressing  will start the stirrer.

The  will flash until the reading is stable. Once the reading is stable you can log and print the measurement by pressing . If you are recording the data into a notebook, record the DO as well as the temperature at which the DO reading is taken. If the stirrer is in motion, press  to turn off the stirrer before removing the stirrer from the sample.

- b. If you are in AUTO-READ™ mode then press  to start the measurement. Once the reading is stable the meter will automatically log and print the results and freeze the display. If the stirrer is enabled, then the stirrer will turn on when  is pressed and turn off once the reading has stabilized.

When using the AUTO-STIR™ BOD probe, pressing the button on the probe will start the AUTO-READ measurement.

- c. If you are in timed measurement mode then the meter will start making measurements as soon as it goes into measurement mode at the frequency selected in setup. It will automatically log and print each measurement. If you are using the benchtop meter and the stirrer control is enabled then pressing  will start the stirrer. Pressing  again will turn off the stirrer.
3. Remove the DO probe from the sample and rinse with deionized water, then place it in your next sample and repeat step 2.
4. When all samples have been measured, rinse the DO probe with deionized water and blot dry. Consult your DO probe manual for proper electrode storage instructions.

Chapter VIII

Conductivity Technique

Conductivity Setup


Description	Display	Range of Values	Default Value	Method Specific
Conductivity Temperature Compensation Selection	C0nd tC LIn	OFF, LIn, nLF Selects temp comp off, linear, or non-linear for natural or ultra-pure water	LIn	Yes
Conductivity Linear Temperature Compensation Coefficient Setting	C0nd COEF 2.1	0.0 - 10.0 Linear temp compensation coefficient in %/C	2.1	Yes
Conductivity TDS Factor Setting	C0nd tdSF 0.49	0.00 - 10.0 TDS factor	0.49	Yes
Conductivity Auto-Calibration Default Cell Constant Setting	C0nd CELL 0.475	0.001 - 199.0 Cell constant used for conductivity AutoCal	0.475	Yes
Conductivity Temperature Reference Selection	C0nd trEF 25	15, 20, 25 Conductivity reference temperature in °C	25	Yes
Conductivity Cell Type & Manual Ranging Selection	C0nd tYPE Std	PLnr, Std, 1, 2, 3, 4, 5, 6, 7 Standard, planar cond cell, or sets manual range 1-7	Std	Yes

1. In the measurement mode press .
2. Press  /  to scroll through the setup menu until **C0nd** is displayed on the top line.
3. Press  to accept the selection and move the arrow icon to the middle line.

4. Press  /  to scroll through the following conductivity settings:
 - **TC** = Temperature Compensation Selection. The users may elect to turn temperature compensation **OFF**, select linear **L In** temperature compensation, or a non-linear **nLF** temperature compensation for natural and ultra-pure water.
 - **COEF** = Temperature coefficient for **L In** (linear) temp compensation expressed in % /C.
 - **TDSF** = TDS Factor for Total Dissolved Solids measurement
 - **CELL** = Conductivity Cell Constant (nominal cell constant used for Auto-Calibration)
 - **TEMP** = Reference temperature for temperature compensation. Options are **15** degrees C, **20** degrees C, or **25** degrees C.
 - **TYPE** = Conductivity cell type (Planar or Conventional)
5. Press  to accept the selection and move the arrow icon to the bottom line.
6. Press  /  and  to edit value.
7. Press  to accept the selection and move the arrow to the top line.
8. Repeat steps 3-7 to edit conductivity settings as needed.
9. Press  to return to measurement mode.

Conductivity Calibration

1. Prepare the probe for use according to the probe user's guide or operator's manual.
2. In the measurement mode press  until the arrow icon is pointing to the conductivity measurement line.
3. Press .
4. Rinse the probe and place into the conductivity standard.
5. To perform a manual calibration – The screen will display the cell constant on the bottom line and the conductivity value on the middle line. Use  and  to change the cell constant until the displayed conductivity value matches the conductivity of the conductivity standard at the measured temperature.
 - If you do not start changing the cell constant within 5 seconds the meter will automatically proceed to the AutoCal™/DirectCal™ screen.
6. To perform an AutoCal or DirectCal calibration – Wait for the **µS/cm** or **mS/cm** icon to stop flashing.
 - AutoCal – When the **µS/cm** or **mS/cm** icon stops flashing the meter will display the temperature corrected value of the conductivity standard.
 - DirectCal – When the **µS/cm** or **mS/cm** icon stops flashing the meter will display the actual conductivity value read by the conductivity probe. Use  and  to change the conductivity value to the actual value of the conductivity standard at the measured temperature.

7. Press  to proceed to the next calibration point and repeat steps 4 through 6 or press  to save and end the calibration.
8. After the value for the last standard has been entered the cell constant will be displayed  appears in the lower field and the actual  constant appears in the main field. The meter will automatically advance to the measurement mode. The  is displayed above the main field.

Conductivity Measurement

1. Rinse conductivity probe and place into sample. Record conductivity directly from the main meter display when the **mS/cm** or **µS/cm** icon stops flashing indicating a stable value. Temperature is displayed in the top left corner of the display.

- 2a. If you are in continuous measurement mode, the instrument will start measuring immediately. If you are using the benchtop meter and the stirrer control is enabled then pressing  will start the stirrer.

The  will flash until the reading is stable. Once the reading is stable you can log and print the measurement by pressing the . If you are recording the data into a notebook, record the conductivity as well as the temperature at which the conductivity reading is taken. If the stirrer is in motion, press  to turn off the stirrer before removing the stirrer from the sample.

- 2b. If you are in AUTO-READ™ mode then press  to start the measurement. Once the reading is stable the meter will automatically log and print the results and freeze the display. If the stirrer is enabled, then the stirrer will turn on when  is pressed and turn off once the reading has stabilized.

- 2c. If you are in timed measurement mode then the meter will start making measurements as soon as it goes into measurement mode at the frequency selected in setup. It will automatically log and print each measurement. If you are using the benchtop meter and the stirrer control is enabled then pressing  will start the stirrer. Pressing  again will turn off the stirrer.

Conductivity Technique

3. Remove the conductivity probe from the sample and rinse with deionized water, then place it in your next sample and repeat step 2.
4. When all samples have been measured, rinse the conductivity probe with deionized water and blot dry. Consult your conductivity probe manual for proper storage instructions.

Chapter IX

ISE Technique

ISE Setup

Description	Display	Range of Values	Default Value	Method Specific
ISE Resolution	ISE rES 1	1, 2, 3 ISE resolution in significant digits	1	Yes
ISE Units	ISE Unlt PPb	M, mG/L, PEr, PPb, nOnE ISE units displayed	PPb	Yes
ISE Calibration Standard Concentration Range	ISE rAng HlgH	LOw, HlgH Stability criteria used during ISE calibration	HlgH	Yes
ISE Auto-Blank Correction	ISE nLI AUtO	AUtO, OFF Enables or disables calibration auto-blank	AUtO	Yes

1. In the measurement mode press .
2. Press  /  to scroll through the setup menu until **ISE** is displayed on the top line.
3. Press  to accept the selection and move to the middle line.

4. Press  /  to scroll through the following options:

- **rES** = Resolution
- **Un It** = ISE Measurement units
 - **M** = Molar
 - **mG/L** = mg/L
 - **PEr** = Percent
 - **PPb** = Parts per billion
 - **nOnE** = No units
- **rAng** = ISE calibration range
 - **HlGH** - for most measurements
 - **Low** - for low level measurements that need more stabilization time
- **nL In** = Non-linear Blank Correction
 - **OFF**
 - **AUtO**

5. Press  to accept the selection and move the arrow icon to the bottom line.

6. Press  /  and  to edit the selected value.

7. Press  to accept the selection and return the arrow icon to the top line.

8. Repeat steps 3 through 7 to edit ISE settings as needed.

9. Press  to return to measurement mode.

Preparation of Standards

The standards should be prepared in the same ISE units as required by the sample results. It is preferable to use serial dilutions with volumetric glassware to obtain the different concentration levels.

Note: *Any reagents such as ionic strength adjustors, should be added to samples and standards as specified in the electrode user's guide or instruction manual. ▲*

The calibration points should bracket the expected concentration range of the samples and there should be a tenfold change in concentration (i.e. 1 ppm and 10 ppm or 10 ppm and 100 ppm).

Fresh aliquots of standard should be used at each calibration.

ISE Calibration

1. Prepare the electrode and standards for use according to the electrode user's guide or instruction manual.
2. Press  until the arrow icon is pointing to the **ISE** measurement line.
3. Press .
4. Rinse the electrode, shake any excess deionized water off and blot dry, and place into the least concentrated standard.
5. Wait for **ISE** to stop flashing. Press  /  and  to change the value of the standard.
6. Press  to proceed to the next lowest calibration standard and repeat steps 4-5, working from lowest concentration to highest, or press  to save the calibration.
7. The slope will be displayed before the meter returns to the measurement mode. **SLP** is displayed in the lower field and the actual electrode slope, in **mV**, is displayed in the main field.

Note: For a 2 or more point calibration the meter will automatically proceed to measurement mode after 3 seconds. For a 1 point calibration the meter will allow you to edit the slope then proceed to measurement mode. ▲

Note: If editing a negative number for slope: ▲

- a. Press  until no digit is blinking and the arrow icon is blinking.
- b. Press  /  to change the sign of the slope.
- c. Press  to save.

ISE Measurement

After calibrating the electrode the meter is ready to start taking measurements.

Make sure that the instrument is in measurement mode ( is lit) and that the **ISE** is lit. If the  is not lit press  to return to measurement mode. If the **ISE** is not lit, press  until the arrow on the left of the display is pointing at the top line. Then press  until the **ISE** lights. You are now ready to start measuring.

1. Rinse the electrode in deionized water. Shake off any excess water and blot dry with lint-free tissue.
2. Place the electrode in your sample.
 - a. If you are in continuous measurement mode, the instrument will start measuring immediately. If you are using the benchtop meter and the stirrer control is disabled, then pressing  will start the stirrer. The **ISE** will flash until the reading is stable. Once the reading is stable you can log and print the measurement by pressing . If the stirrer is on, press  to turn off the stirrer before removing the stirrer from the sample.
 - b. If you are in AUTO-READ™ mode then press  to start the measurement. Once the reading is stable the meter will automatically log and print the results and freeze the display. If the stirrer is enabled then the stirrer will turn on when  is pressed and turn off once the reading has stabilized.
 - c. If you are in timed measurement mode then the meter will start making measurements at the frequency you have specified during setup. It will automatically log and print each measurement. If you are using the benchtop meter and the stirrer control is enabled, then pressing  will start the stirrer. Pressing  again will turn off the stirrer.

ISE Technique

3. Remove the electrode from the sample and rinse with deionized water, shake or blot dry, then place it in your next sample and repeat step 2.
4. Once you are finished measuring samples, rinse the electrode with deionized water and blot dry. Then consult your electrode manual for proper electrode storage instructions.

Chapter X

Terms and Conditions

General

Seller hereby offers for sale to the buyer ("Buyer") the products herein (the "Products") on the express condition that Buyer agrees to accept and be bound by the terms and conditions set forth herein. Any provisions contained in any document issued by Buyer are expressly rejected and if the terms and conditions in this Agreement differ from the terms of Buyer's offer, this document shall be construed as a counter offer and shall not be effective as an acceptance of Buyer's document. Buyer's receipt of Products or Seller's commencement of the services provided hereunder will constitute Buyer's acceptance of this Agreement. This is the complete and exclusive statement of the contract between Seller and Buyer with respect to Buyer's purchase of the Products. No waiver, consent, modification, amendment or change of the terms contained herein shall be binding unless in writing and signed by Seller and Buyer. Seller's failure to object to terms contained in any subsequent communication from Buyer will not be a waiver or modification of the terms set forth herein. All orders are subject to acceptance in writing by an authorized representative of Seller.

Warranty

The warranty for Thermo Scientific Orion products covers failures due to manufacturer's workmanship or material defects from the date of purchase by the user. User should return the warranty card and retain proof of purchase. Warranty is void if product has been abused, misused, or repairs attempted by unauthorized persons.

Warranties herein are for product sold/installed by Thermo Fisher Scientific or its authorized dealers.

Any product sold by a U.S. or Canadian distributor must be returned to Thermo Fisher Scientific for any warranty work. Please contact our Technical Service department for further information. A Return Authorization Number must be obtained from Technical Service before returning any product for in-warranty repair

or replacement. In the event of failure within the warranty period Thermo Fisher Scientific will, at Thermo Fisher Scientific's option, repair or replace product not conforming to this warranty. There may be additional charges, including freight, for warranty service performed in some countries. For service, call Thermo Fisher Scientific (or its authorized dealer outside the United States and Canada). Thermo Fisher Scientific reserves the right to ask for proof of purchase, such as the original invoice or packing slip.

Field Service is available on the BOD AutoEZ™. Contact our Field Service department for details on quotations, service, and other field service-related activities.

The following products are warranted to be free from defects in material and workmanship in the period listed below from the date of purchase from the user or from the date of shipment from Thermo Fisher Scientific, whichever is earlier, provided use is in accordance with the operating limitations and maintenance procedures in the instruction manual and when not having been subjected to accident, alteration, misuse, abuse or breakage of electrodes:

Thirty-six months from date of purchase by the user (or forty-two months from date of shipment from Thermo Fisher Scientific)

All Orion Star™ Series pH, ISE, DO and Conductivity Meters, Waterproof Meters (630, 635, 830A, 835A, 260A, 261S, 265A, 266S, 130A, 131S, 135A, 136S), Conductivity Meters (105Aplus™, 115Aplus™, 125Aplus™, 145Aplus™, 150Aplus™ and 162A), PerpHect® pH and pH/ISE Meters (310, 320, 330, 350, 370), pH and pH/ISE Meters (210Aplus™, 230Aplus™, 250Aplus™, 290Aplus™, 410Aplus™, 420Aplus™, 520Aplus™, 525Aplus™, 710Aplus™, 720Aplus™ and 920Aplus™), pHuture MMS™ Meters (535A and 555A), pH/Conductivity Meter (550A), Dissolved Oxygen Meters (805Aplus™, 810Aplus™, 850Aplus™ and 862A).

Twenty-four months from date of purchase by the user (or thirty-six months from date of shipment from Thermo Fisher Scientific)

Ross Ultra® Electrodes, AQUAfast® IV Colorimeters, AQUAfast® IV Turbidimeter, 925 Flash Titrator™, Series 100 DuraProbe™ Conductivity Cells and Series 800 Dissolved Oxygen Probes.

Twelve months from date of purchase by the user (or eighteen months from date of shipment from Thermo Fisher Scientific)

Laboratory pH Meters, (301, 611 and 940), SensorLink®, pHuture® pH Meters (610 and 620), Smart Chek™ Meters, Sage® Pumps, Cahn® Balances, 930 Ionalyzer®, 950 ROSS® FAST QC™ Titrator, 960 Titrator PLUS®, Karl Fischer Titrators, Autosamplers, BOD AutoEZ™, pHuture® Conversion Box, Wine Master®, 607 Switchbox, rf link™, AQUAfast® II Colorimeters.

Ion Selective Electrodes, ionplus® Electrodes, ROSS® Electrodes, Sure-Flow® Electrodes, PerpHecT® Electrodes, AquaPro Professional Electrodes, No Cal® pH Electrodes, Standard pH Electrodes, TRIS pH Electrodes, KNipHE® Electrodes, ORP Triode™ (9180BN), pHuture® pH Probe (616500) and pHuture MMS™ Quatrode™ and Triode™ (616600 and 617900), 9708 DO Probe, Series 100 Conventional Conductivity Cells, Temperature Probes and Compensators (except those models noted).

The 93 and 97 ionplus® Series Sensing Modules are warranted to give six months of operation if placed in service before the date indicated on the package, except 9307 and 9707 Nitrate Modules are warranted to give ninety days of operation if placed in service before the date indicated on the package.

Six months from date of purchase by the user (or twelve months from date of shipment from Thermo Fisher Scientific)

Flash Titration™ Probe (092518), pHuture® Electrode (615700), pHuture MMS™ Pentrode™ (617500), Quatrode™ (617800) and Triode™ (615800), Low Maintenance Triode™ (9107BN), ORP Low Maintenance Triode™ (9179BN), and PerpHecT® Low Maintenance Triode™ (9207BN), Waterproof Triode™ (9107WP, 9107WL, 9109WL and 9109WP), QuiKcheK® Meters and Micro Electrodes.

Three months from date of purchase by the user (or six months from date of shipment from Thermo Fisher Scientific)

Economy Line Electrodes, 9105, 9106, 9115, 9116, 9125, 9126, 9135, 9136, 9206. Warranty also includes failure for any reason (excluding breakage), except abuse, provided the electrode is not used in solutions containing TRIS, silver, sulfide, perchlorate, or hydrofluoric acid; or in solutions more than one (1) Molar in strong acid or base at temperatures above 50 °C.

“Out-of-Box” Warranty - Should any of the following products fail to work when first used, contact Thermo Fisher Scientific immediately for replacement.

Solutions, Standards, Reagents, Cables, Line adapters, Printers, Software, Cases, Stands, Probe Membranes, AQUAfast® Test Strips, and general accessories.

For products in the catalog not listed in this warranty statement, please visit our website at: www.thermo.com/water.

THE WARRANTIES DESCRIBED ABOVE ARE EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES WHETHER STATUTORY, EXPRESS OR IMPLIED INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE AND ALL WARRANTIES ARISING FROM THE COURSE OF DEALING OR USAGE OF TRADE. THE BUYER'S SOLE AND EXCLUSIVE REMEDY IS FOR REPAIR OR REPLACEMENT OF THE NON-CONFORMING PRODUCT OR PART THEREOF, OR REFUND OF THE PURCHASE PRICE, BUT IN NO EVENT SHALL THERMO FISHER SCIENTIFIC (ITS CONTRACTORS AND SUPPLIERS OF ANY TIER) BE LIABLE TO THE BUYER OR ANY PERSON FOR ANY SPECIAL, INDIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES WHETHER THE CLAIMS ARE BASED IN CONTRACT, IN TORT (INCLUDING NEGLIGENCE), OR OTHERWISE WITH RESPECT TO OR ARISING OUT OF THE PRODUCT FURNISHED HEREUNDER. REPRESENTATION AND WARRANTIES MADE BY ANY PERSON, INCLUDING ITS AUTHORIZED DEALERS, REPRESENTATIVES AND EMPLOYEES OF THERMO FISHER SCIENTIFIC WHICH ALTER OR ARE IN ADDITION TO THE TERMS OF THIS WARRANTY SHALL NOT BE BINDING UPON THERMO FISHER SCIENTIFIC UNLESS IN WRITING AND SIGNED BY ONE OF ITS OFFICERS.

Limitation of Liability

NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, THE LIABILITY OF SELLER UNDER THESE TERMS AND CONDITIONS (WHETHER BY REASON OF BREACH OF CONTRACT, TORT, INDEMNIFICATION, OR OTHERWISE, BUT EXCLUDING LIABILITY OF SELLER FOR BREACH OF WARRANTY (THE SOLE REMEDY FOR WHICH SHALL BE AS PROVIDED UNDER SECTION 2 ABOVE)) SHALL NOT EXCEED AN AMOUNT EQUAL TO THE LESSER OF (A) THE TOTAL PURCHASE PRICE THEREFORE PAID BY BUYER TO SELLER WITH RESPECT TO THE PRODUCT(S) GIVING RISE TO SUCH LIABILITY OR (B) ONE MILLION DOLLARS (\$1,000,000). NOTWITHSTANDING ANYTHING TO THE CONTRARY CONTAINED HEREIN, IN NO EVENT SHALL SELLER BE LIABLE FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES (INCLUDING WITHOUT LIMITATION DAMAGES FOR LOSS OF USE OF FACILITIES OR EQUIPMENT, LOSS OF REVENUE, LOSS OF DATA, LOSS OF PROFITS OR LOSS OF GOODWILL), REGARDLESS OF WHETHER SELLER (a) HAS BEEN INFORMED OF THE POSSIBILITY OF SUCH DAMAGES OR (b) IS NEGLIGENT.

Miscellaneous

(a) The rights and obligations of the parties hereunder shall be governed by and construed in accordance with the laws of the Commonwealth of Massachusetts, without reference to its choice of law provisions. Each party hereby irrevocably consents to the exclusive jurisdiction of the state and federal courts located in Suffolk County, Massachusetts, in any action arising out of or relating to this Agreement and waives any other venue to which it may be entitled by domicile or otherwise. (b) In the event of any legal proceeding between the Seller and Buyer relating to this Agreement, neither party may claim the right to a trial by jury, and both parties waive any right they may have under applicable law or otherwise to a right to a trial by jury. Any action arising under this Agreement must be brought within one (1) year from the date that the cause of action arose.

NOTES

Chapter XI

Declaration of Conformity

Manufacturer: Thermo Fisher Scientific Inc.

Address: 166 Cummings Center
Beverly, MA 01915
USA

We declare that the following products described below conform to the Directive and Standard listed below:

Product(s): Meters for Measuring pH, Conductivity, Dissolved Oxygen, and/or ISE
Benchtop models are rated 100-240 VAC, 50/60 Hz, 0.5 A
Handheld models use 4 non-rechargeable AA Batteries

Benchtop	Portable
5-Star pH/ISE/Conductivity/DO Meter	5-Star pH/ISE/Conductivity/DO Meter
4-Star pH/Conductivity Meter	5-Star pH/Conductivity/DO Meter
4-Star pH/DO	4-Star pH/Conductivity Meter
4-Star pH/ISE Meter	4-Star pH/DO Meter
3-Star Conductivity Meter	4-Star pH/ISE Meter
3-Star DO Meter	3-Star Conductivity Meter
3-Star pH Meter	3-Star pH Meter
2-Star pH Meter	3-Star DO Meter

Equipment Class: Measurement, control and laboratory
Benchtop models are EMC Class A
Portable models are EMC Class D

Declaration of Conformity

Directive(s) and Standard(s):

- 89/336/EEC – Electromagnetic Compatibility (EMC Directive)
 - EN 61326:1997 + A1:1998 + A2:2001 – Electrical equipment for measurement, control, and laboratory use – EMC requirements
- 73/23/EEC – Low Voltage Directive (LVD)
 - EN 61010-1:2001 – Safety requirements for electrical equipment for measurement, control, and laboratory use – general requirements

Manufacturer's Authorized Representative:

Date:

Steve Brody
Manager, Compliance Engineering
Thermo Fisher Scientific Inc.

June 6, 2006

WEEE Compliance:

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC .
It is marked with the following symbol:

Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State and this product should be disposed of or recycled through them. Further information on compliance with these Directives, the recyclers in your country, and information on Thermo Scientific Orion products which may assist the detection of substances subject to the RoHS Directive are available at www.thermo.com/WEEERoHS.

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